

# A Stem Cell Potency and Release Criteria Assay Specifically Designed for Umbilical Cord Blood Transplantation that is Compliant with Regulatory Guidelines


Ivan N. Rich, Karen M. Hall and Holli Harper  
**HemoGenix®**, Inc. Colorado Springs, CO,

## POTENCY IS A MEASURE OF CHEMICAL OR BIOLOGICAL ACTIVITY. HOW IS POTENCY MEASURED?

"All potency assays used for release testing of licensed biological drug products must comply with applicable biologics and cGMP regulations including:

1. Indicate potency (biological activity/activities) specific to the product.
2. Provide test results for release of the product.
3. Meet pre-defined acceptance and/or rejection criteria.
4. Include appropriate reference materials, standards, and/or controls.
5. Establish and document the accuracy, sensitivity, specificity, reliability, reproducibility and robustness of the test methods, employed through validation".

(FDA Draft Guidance for Industry: Potency Tests for Cellular and Gene Therapy Products, October 2008). To measure potency, a reference standard is an absolute requirement. A dose response is performed for the reference standard and the samples. Comparison of the dose response lines should demonstrate statistical parallelism as seen in the example below. The horizontal displacement from the reference standard provides the potency ratio and therefore the activity of the sample.



NO test or assay (TNC, viability, CD34 or CFC) presently used by stem cell transplantation and cord blood storage processing laboratories complies with these criteria (FDA and EMEA) for a potency assay designed for the intended use of the product.

## Umbilical Cord Blood Stem Cell Potency, Release and Engraftment Study

1. Performed in collaboration with the University of Colorado Cord Blood Bank (ClinImmune, Inc).
2. 56 Umbilical cord blood frozen pellet samples for which engraftment information was known.
3. Approx. 700µl / sample.
4. All samples tested randomly.
5. Nycoprep MNC, cell count, viability by 7-AAD.
6. 24 Samples tested for CFC-GEMM and 23 samples tested for HPP-SP at 2,500, 5,000 and 7,500 or 10,000 cells/well, when sufficient cells were available.
7. All remaining samples tested at 5,000 cells/well.

## Characteristics of the Study Population

Time from freezing to transplantation: 91 days to 9 years and 311 days.  
 Time from transplantation to ANC recovery (> 500 cells/µl): 5 - 114 days.  
 Time from transplantation to platelet recovery (> 50k/µl): 2 - 237 days.

## Cord Blood Characteristics


Pre-cryopreservation total nucleated cell count:  $15.7 \times 10^6$  -  $178 \times 10^6$ .  
 Post-thaw total nucleated cell count: 20,000 -  $5.55 \times 10^6$ .  
 Post-thaw viability (7-AAD): 49% - 100%.

## Cord Blood Storage, ATP Values and Engraftment


32 of 56 UCB samples stored at -80°C (but units stored in LN2).  
 17 out of these 32 units engrafted, but ALL 32 samples tested exhibited ATP values <0.04µM.

24 of 56 samples (and units) stored in LN2.  
 ALL 24 engrafted (100%) and ALL exhibited ATP values >0.04µM for CFC-GEMM or >0.05µM for HPP-SP or both.


## Important Correlations (or lack of) Associated with ATP Concentration and Stem Cell Potency



## Neither the CFC-GEMM nor HPP-SP (as concentration or per kilogram transplanted) Correlates with Neutrophil Engraftment (similar results for platelets, no shown)



## Measurement of Cord Blood Stem Cell Potency



Of the 24 UCB samples that had been stored in LN2 and for which a 3-point cell dose response could be performed, only 3 UCB samples produced dose response curves statistically parallel to the UCB reference standard. The common CFC-GEMM slope for UCB 8 and 14 was  $0.97 \times 10^5$ , while that for HPP-SP was  $2.93 \times 10^5$ . The CFC-GEMM potency ratio for UCB 8 was 0.94 and that for UCB 14 was 1.05. The HPP-SP potency ratio for UCB 3 was 1.36. A potency ratio greater than 1, indicates a greater potency than the reference standard.

## What does the correlation of slope with iATP concentration mean?

- The steeper the slope of the cord blood dose response:
1. The greater the proliferation potential.
  2. The more primitive the stem cells in the sample.
  3. The greater the stem cell potency, and
  4. The greater the probability of engraftment.

## REFERENCES

1. Development of a Novel Assay to Evaluate the Functional Potential of Umbilical Cord Blood Progenitors, Reems, Hall, Ludalay, Taber, Rich, Transfusion, 48:620-628 (2008).
2. Cell Potency Assays for the 21st Century Stem Cell Transplantation and Cord Blood Bank Processing Laboratories. HemoGenix® White Paper. Please download this article from the HemoGenix® website at [www.hemogenix.com](http://www.hemogenix.com).
3. "Validation and Development of a Predictive Paradigm for Hemotoxicity Using a Multifunctional Bioluminescence Colony-Forming Proliferation Assay". Rich IN, Hall KM. Tox Sci 87: 427-441 (2005).
4. In vitro to in vivo concordance of a high throughput assay of bone toxicity across a diverse set of drug candidates. Olaharski et.al. Tox. Lett. 2009. Science Direct Online.

## HALO®-96 PQR

An instrument-based, ATP bioluminescence proliferation assay that measures cord blood stem cell potency and defines release criteria.

**ATP - A Powerful Biochemical Marker of Cellular Function**

**The HALO® Principle**

As cells proliferate, intracellular ATP increases. The HALO®-96 PQR measures intracellular ATP (iATP) using a proprietary ATPase enzyme that converts ATP to ADP and releases Pi. The released Pi is then converted to ATP by a luciferase enzyme, which produces light. The intensity of the light is proportional to the amount of iATP present.

HALO®-96 PQR utilizes a proprietary Suspension Expansion Culture (SEC) technology, rather than methyl cellulose. This allows greater ease of use, rapid 5 day assay completion, increased sensitivity and low coefficients of variation (see below). Unlike the CFC assay, HALO®-96 PQR (like all HALO® products) is standardized and validated according to FDA regulations and guidelines. All HALO® products show equivalency to the CFC assay and therefore can replace the CFC assay for most applications.

**HALO®-96 PQR Precision: Part of the validation process to test reproducibility**

Background Controls				
Cells / Well	1,000	1,000	1,000	10,000
No. of Cultures	10	10	10	10
Mean CV	5.8	6.3	3.7	3.4


CFC-GEMM				
Cells / Well	1,000	1,000	1,000	10,000
No. of Cultures	10	10	10	10
Mean CV	16.6%	16.4%	9.4%	8.6%

HPP-SP				
Cells / Well	1,000	1,000	1,000	10,000
No. of Cultures	10	10	10	10
Mean CV	12.8%	9.7%	6.7%	5.7%

**Correlations:**  
 Cell dose vs CFC: R = 0.99  
 Cell dose vs HALO®-96 SEC: R = 0.987  
 HALO®-96 SEC vs CFC: R = 0.986.

**Conclusion:**  
 HALO®-96 can replace the CFC Assay

## Relationship between iATP, Release Criteria, Stem Cell Potency and Engraftment



## SUMMARY

1. TNC, viability, CD34 and CFC cannot be used as potency assays because they are not compliant with the regulations for potency.
2. iATP concentrations can distinguish between proliferative and non-proliferative cells.
3. Cellular potency must take the biology and physiology of the cells into account.

## CONCLUSIONS

1. HALO®-96 PQR is a rapid, reference standard-based stem cell potency assay for umbilical cord blood that can help define acceptance limits for release criteria.
2. HALO®-96 PQR is validated and fully compliant with FDA and EMEA regulations and guidelines.
3. Stem cell potency is determined by the cumulative stem cell proliferation potential.